

**POSITION PERFECT™**  
POSITIONING AIDS




**CONTOURED FOR PATIENTS  
AND TECHNOLOGISTS™**

**contour fabricators inc.**

## **2800 Urology Table Pad**

Contour Fabricators, Inc. introduces the 2800 Urology Table Pad. Designed specifically for use with the GE/OEC UroView 2800, the pad's unique scooped design and attached flap direct fluids away from the patient and into the drain bag. The durable Scancoat™ provides a cleanable, stain resistant barrier.



**Compatibility:**

Urology

**Applications:**

Designed Specifically for GE/OEC UroView 2800

**Dimensions:**

30.5" x 26" x 2.5"

**Product Features:**

Flap

**Benefits:**

Directs fluids into drain bag

Contoured Groove

Directs fluids away from patient to the base of the pad

**Related Products**

CFI custom builds Urology Table Pads to fit any Urology Table

Phone (810) 750-5300  
Fax (810) 750-5310

**Contour Fabricators**  
14241 Fenton Road  
Fenton, MI 48430

sales@contourfab.com  
www.contourfab.com