



## contour fabricators inc.

# Camera Drapes

Contour Fabricators, Inc. introduces its line of custom Camera Drapes. Used to keep the camera sterile during surgical procedures, the drapes are made out of polyethylene and are telescope folded with a cuff to allow for quick insertion. All have double sided tape for the gathering of material and secure closure. They are latex free with easy to follow instruction labels.



**Compatibility:**

All camera types

**Dimensions:**

Part # 5431 - 10" x 80", clear polyethylene with perforated tip

Part # 5432 - 5" x 96", clear polyethylene with perforated tip

Part # 5433 - 7" x 96", clear polyethylene with elasticized fenestration to provide a secure fit around the hand piece

Part # 1500 - 5" x 96", blue camera cover with rubber band to provide a secure fit around the hand piece

**Options:**

Camera Drapes can be customized to fit any camera

**Related Products**

CFI manufactures other drapes and covers including Microscope Drapes, C-Arm Drapes, Band Bags and Dome Bags