


E.I.B.O.S. Protective Drape

Contour Fabricators, Inc. presents the E.I.B.O.S. Protective Drape. Designed to keep the Moeller E.I.B.O.S. attachment sterile during surgical procedures, the E.I.B.O.S. Protective Drape's contoured design allows it to fit easily and securely into position around the E.I.B.O.S. attachment.



The drape is made of clear polyethylene with a black polyethylene frame. The E.I.B.O.S. protective drape includes a foam ring used for integration with the full microscope drape, a polycarbonate lens, and a PVC end cap.

Compatibility:

Moeller E.I.B.O.S Microscope

Applications:

All E.I.B.O.S microscope procedures

Dimensions:

7" x 12"

Product Features:

Benefits:

Contoured Fit

Reduces excess material around E.I.B.O.S attachment

Foam Ring

Allows for 360 degree rotation and a 90 degree pivot
Provides easy integration with full microscope drape

Clear Polyethylene

Creates a completely sterile field around microscope

Polycarbonate Lenses

Provide enhanced lens clarity

Sterilizable bottom lens end cap

Contributes to sterility of microscope