



contour fabricators inc.

Femoral Weights

Contour Fabricators, Inc. presents the Femoral Weights. Designed to maintain pressure on the patient following vascular/catherization studies, the Femoral Weights provide enough pressure to eliminate the need for clamping devices or pressure application by hospital staff.



Compatibility:

Contain lead - not compatible in imaging.

Applications:

Catheterization Studies

Dimensions:

CFI-505 – 3" x 9", 5 lb.

Part # 21934 – 5" x 9", 8 lb.

Part # 20659 – 6" x 17", 10 lb.

Product Features:

Heavy and narrow

Durable, double layered vinyl

Benefits:

Provides localized pressure and contours to body

Will not leak, easily cleaned and easily taped in place

Options:

CFI manufactures a complete line of weighted immobilization and compression devices in a variety of weights and materials.