



contour fabricators inc.

Microscope Drapes

Contour Fabricators introduces its line of custom Microscope Drapes designed for the Zeiss, Moeller and Leica Microscopes.



The drapes are made out of a strong, drapable polyethylene film to prevent tearing and provide a completely sterile field. They also include monocular or stereo observer tubes for enhanced vision and straps to gather the excess poly and the re-attachable fiber-free ocular closures. The drapes can be custom fit to the lengths of any Microscope.



Compatibility:

Zeiss, Moeller and Leica Microscopes

Dimensions:

Zeiss MD & Moeller

part # 441265 - 44" wide x 120" long, fits 65 mm lens

part # 421265 - 42" wide x 120" long, fits 65 mm lens

part # 446565 - 44" wide x 65" long, fits 65 mm lens

Zeiss OPMI

part # 446050 - 44" wide x 150" long, fits 50 mm lens

Leica Wild

part # 541568 - 54" wide x 150" long, fits 68 mm lens

Options:

Microscope Drapes can be customized to fit any microscope

Related Products

CFI manufactures other drapes and covers including Camera Drapes, C-Arm Drapes, Band Bags and Dome Bags