



contour fabricators inc.

Performance Plus[®]



c-arm covers • microscope and camera drapes • miscellaneous disposable covers and bags • urology disposables • foam positioning aids

The
Power
of **Selection**

Mini C-Arm Covers

Part #	Description of Part	Unit of Measure	Detailed Description
5400	Mini C-Arm Kit, fits OEC Mini 6600 with grey footswitch	20/case	Kit consists of a Mini C-Arm Drape, Mini-C Keyboard Cover, and a 17" x 15" Footswitch Cover. Sterile, individually pouched.
5401	Mini C-Arm Kit, fits OEC Mini 6600 with black footswitch	20/case	Kit consists of a Mini C-Arm Drape, Mini-C Keyboard Cover and a 10" x 8" Footswitch Cover. Sterile, individually pouched.
5402	Mini C-Arm Kit, fits OEC Mini 6800	20/case	Kit consists of a MiniView C-Arm Drape, Keyboard Cover and a 17" x 15" grey Footswitch Cover. Sterile, individually pouched.
5403	Mini C-Arm Drape, fits OEC Mini 6800	25/case	Single piece drape to cover I.I., X-ray Tube and the C-Arm. Sterile, individually pouched.
5404	Mini C-Arm Drape, fits OEC Mini 6600	20/case	Single piece drape to cover I.I., X-ray Tube and the C-Arm. Sterile, individually pouched.
5423	Mini C-Arm Drape, Fluoroscan Premier	20/case	Designed as one piece to fit the I.I., X-ray Tube and C-Arm. Sterile, individually pouched.
5405	Mini C Keyboard Cover	25/case	Clear poly, 14" x 17" with double sided tape. Sterile, individually pouched.



C-Arm Covers

Part #	Description of Part	Unit of Measure	Detailed Description
5428	9" I.I. C-Arm Kit, Fits OEC 9600/9800	20/case	Consists of C-Arm Drape with 5 mounting clips, 30" x 32" I.I. Cover, 30" x 32" X-Ray Tube Cover, and a 17" x 15" Footswitch Cover. Sterile.
5429	12" I.I. C-Arm Kit, Fits OEC 9600/9800	20/case	Consists of C-Arm Drape with 5 mounting clips, 36" x 36" I.I. Cover, 36" x 36" X-Ray Tube Cover and a 17" x 15" Footswitch Cover. Sterile.
5427	C-Drape Fits OEC 9600/9800	30/case	C-Arm Drape with 5 mounting clips. Sterile, individually pouched.
5411	9" I.I. / X-Ray Tube Cover	60/case	30" x 32" Bandbag. Sterile, individually pouched.
5412	12" I.I. / X-Ray Tube Cover	30/case	36" x 36" Bandbag. Sterile, individually pouched.
5413	IR Hand Control Cover	25/case	7" x 13" clear cover. Sterile, individually pouched.
5414	9" I.I. Universal Halfbag Cover	20/case	36" x 80" designed to fit the I.I. and midway down the C-Arm. Sterile, individually pouched.



Miscellaneous Disposable Covers


Part #	Description of Part	Unit of Measure	Detailed Description
CFI-967	Mayo Stand Cover 20" X 54"	25/case	Sterile, individually pouched.
CFI-708	Mayo Stand Cover 24" X 58"	25/case	Sterile, individually pouched.
CFI-966	Mayo Stand Cover 36" X 72"	25/case	Sterile, individually pouched.
CFI-939	Set-Up Cover 36" X 84"	100/case	With caution label. Sterile, individually pouched.
CFI-702	Set-Up Cover 44" X 70"	25/case	With caution label. Sterile, individually pouched.
CFI-874	Set-Up Cover 54" X 70"	100/case	With caution label. Sterile, individually pouched.
CFI-704	Slush Drape 54" X 54"	25/case	Sterile, individually pouched.
CFI-1501	Blue Ultrasound Probe Cover 5" X 72"	20/case	Sterile, individually pouched.



Microscope and Camera Drapes


Drapes for Zeiss MD® Microscope

Part #	Description of Part	Unit of Measure	Detailed Description
441265	44.0" wide x 120.0" long Fits 65mm lens	10/case	All of our microscope drapes are made with strong, drapable polyethylene film, monocular or stereo observer tubes, and include straps to gather excess poly and re-attachable fiber-free ocular closures.
421265	42.0" wide x 120.0" long Fits 65mm lens	10/case	
446565	44.0" wide x 65.0" long Fits 65mm lens	10/case	




Drapes for Zeiss OPMI® Microscope

446050	44.0" wide x 60.0" long Fits 50mm lens	10/case	Contoured to reduce unnecessary excess drape.
--------	---	---------	---




Drapes for Leica Wild® Microscope

541568	54.0" wide x 150.0" long Fits 68mm lens	10/case	
--------	--	---------	--



Camera Drapes

5431	Camera Drape, 10" x 80", clear polyethylene with perforated tip.	25/case	Camera drapes are telescope folded with cuff to allow for quick insertion. All have double-sided tape for gathering of material and secure closure. Latex-free with easy to follow instruction labels. Sterile, individually pouched.
5432	Camera Drape, 5" x 96", clear polyethylene with perforated tip.	25/case	
5433	Camera Drape, 7" x 96", clear polyethylene with elasticized fenestration to provide a secure fit around the hand piece.	25/case	
CFI-1500	Blue Camera Cover 5" X 96"	20/case	



Miscellaneous Disposables

Band Bags – Clear poly, with a latex-free rubber band heat sealed in one end. Sterile, individually pouched.

Part #	Part Dimensions	Unit of Measure
CFI-654	36" X 20"	100/case
CFI-655	36" X 30"	100/case
CFI-657	36" X 36"	100/case
CFI-667	30" X 36"	100/case
CFI-761	40" X 30"	100/case
CFI-827	40" X 40"	100/case



Dome Bags – Shower cap style bag, latex-free elastic is sewn around the perimeter of the material. Sterile, individually pouched.

Part #	Part Dimensions	Unit of Measure
CFI-756	27" depth	50/case
CFI-821	22" depth	50/case
CFI-893	26" depth	50/case
CFI-911	18" depth	100/case



Foot Switch Covers – Clear poly, non-sterile, different sizes, latex-free elastic cord with cinch lock.

Part #	Part Dimensions	Unit of Measure
CFI-706	12" X 12"	50/case
CFI-962	20" X 12"	50/case
CFI-707	12" X 8"	50/case
CFI-963	24" X 20"	50/case
CFI-961	14" X 30"	50/case
CFI-964	30" X 30"	50/case
5406	17" x 15"	100/case
5407	10" x 8"	100/case
5421	24" x 20"	100/case
5422	20" x 36"	100/case



Used with OEC Grey Footswitch.

Wide Footswitch Cover fits with cinch string closer. OEC UroView 2600 Wide Footswitch bulk packaged. Non-sterile.

Wide Footswitch Cover fits OEC with cinch string closer. UroView 2800 Wide Footswitch Bulk packaged. Non-sterile.

Cassette Covers Sterile and Non-Sterile – Clear poly sterile covers are individually pouched. Non-sterile are bulk packed 100 per case.

Part #	Part Dimensions	Unit of Measure
CFI-954	8" X 10"	25/case
CFI-955	10" X 12"	25/case
CFI-956	14" X 17"	25/case
CFI-957	18" X 20"	25/case
CFI-958	21" X 40"	25/case
CFI-959	24" X 24"	25/case
CFI-705	24" X 36"	25/case
CFI-960	30" X 50"	25/case

Sterile
Sterile
Sterile
Sterile
Sterile
Sterile
Sterile
Sterile



Drain Bag Frames and Accessory Brackets

Part #	Description of Part	Unit of Measure	Detailed Description
21848	Urology Hoop Frame	each	Mounting brackets and hoop frame are sold as a set. Designed to fit all surgical tables including Steris, Skytron, and Amsco. Quick and easy installation.
22672	Urology Frame and Brackets Kit	each	Bracket and hoop frame integrated as one piece. Quick and easy installation.



Urology Disposables

Part #	Description of Part	Unit of Measure	Detailed Description
5313	Uro-Catcher Bag, Fits all hoop style frames	20/case	Designed to fit all hoop frames. Vinyl flap to go under the patient to direct fluids into the drain bag. Tissue filter and a flexible hose that extends 6 feet. Sterile, individually pouched.
5415	Urology Disposables Pack, Fits OEC UroView 2600	20/case	Pack consists of 5419 Drain Bag, Hand Control Cover, 17" x 15" Footswitch Cover and a 24" x 20" Footswitch Cover. Sterile, individually pouched.
5416	Urology Disposables Pack, Fits OEC UroView 2800	20/case	Pack consists of 5419 Drain Bag, Hand Control Cover, 17" x 15" Footswitch Cover and a 20" x 36" Footswitch Cover. Sterile, individually pouched.
5417	Square Drain Bag with plastic tabs, injection molded screen, and funnel. Fits all square frames.	20/case	Both bags 5417 and 5000 fit the GE/OEC Uroview 2000 and GE/OEC Uroview 2500, Liebel Flarsheim/ Mallinckrodt Hydrovision ED, DR, and HP Urology Systems with the collapsible bracket, fixed bracket, or square bracket. Sterile, individually pouched.
5000	Square Drain Bag with metal tabs, and mesh screen. Fits all square frames.	20/case	
5419	Hoop Drain Bag, Fits OEC 2600/2800	20/case	Designed to fit the GE/OEC Uroview 2600 and GE/OEC Uroview 2800 tables. Heavy gauge vinyl strip to hold bag snugly to the foot of the table. Mesh screen and a flexible hose that extends 6 feet. Sterile, individually pouched.
5010	Urology Drain Collection Bag	10/case	Designed to collect waste fluid when a floor drain or suction is not available. Internal baffle system enables the container to fill evenly as powder converts liquids to gel. Built in handles to transport the unit for disposal. 5 gallon capacity. Patent # 5,549,707
5420	Hand Control Cover	100/case	5" x 10" with rubber band closer. Bulk packaged. Non-sterile



Urology Table Pads

Part #	Description of Part	Unit of Measure	Detailed Description
21443	2500/2600 Soft Table Pad Low density open cell foam	each	Scooped design with flap to help guide and channel fluid into drain bag, coated.
22423	2500/2600 Firm Table Pad High density open cell foam	each	
22421	2500/2600 Soft Table Extension Pad	each	Low density open cell foam, coated.
22422	2500/2600 Firm Foot Extension Pad	each	High density open cell foam, coated.
22710	2800 Standard Table Pad	each	Medium density open cell foam, coated. Scooped design with flap to help guide and channel fluids into the drain bag.
22708	2800 Premium Table Pad	each	Closed cell foam, coated. Scoop design with flap to help guide and channel fluids into the drain bag.
22711	2800 Standard Foot Extension Pad	each	Medium density open cell foam, coated.
22709	2800 Premium Foot Extension Pad	each	Closed cell foam, coated.
20983	Urology Table Pad Liebel Flarsheim	each	Medium density open cell foam, coated. Scooped design with flap to help guide and channel fluids into the drain bag.
22692	Foot Extension 16" Liebel Flarsheim	each	Medium density open cell foam, coated.
22691	Foot Extension 24" Liebel Flarsheim	each	Medium density open cell foam, coated.



Procedural Foam - Disposable Foam Positioning Aids

Part #	Description of Part	Unit of Measure	Detailed Description
CFI-644	Circular Head Rest, 2" x 7" diameter	25/case	Blue foam, bulk packed, non-sterile.
CFI-645	Circular Head Rest, 2" x 9" diameter	20/case	Blue foam, bulk packed, non-sterile.
CFI-680	Sculpted Head Rest, 7" x 7" x 2"	20/case	Blue foam, bulk packed, non-sterile.
CFI-866	Concave Disc, 1.75" x 7" diameter	20/case	Blue foam, bulk packed, non-sterile.
CFI-867	Concave Disc, 2.50" x 7" diameter	20/case	Blue foam, bulk packed, non-sterile.
CFI-868	Concave Disc, 3.00" x 7" diameter	20/case	Blue foam, bulk packed, non-sterile.
CFI-869	Cervical Head Rest 3.5" x 4.5" x 8"	each	Blue foam, non-sterile.
CFI-870	Spinal Wedge 7" x 7" x 21.25"	each	Blue foam, non-sterile.
CFI-871	Pressure Block 5.25" x 7" x 10"	each	Blue foam, non-sterile.
CFI-734	Cylinder 5" diameter x 8"	20/case	Blue foam, bulk packed, non-sterile.
CFI-875	Cylinder 5" diameter x 16"	10/case	Blue foam, bulk packed, non-sterile.

