


Secure Slide System™

Patent Pending

Contour Fabricators, Inc. presents the Secure Slide System™. No longer are technicians and practitioners required to struggle with positioning and securing voluntary patient immobilization straps. Our patented *slide groove* approach allows the straps and other attachments to slide freely along the length of the table, providing infinite adjustability. The convenience of this system can be adapted to any imaging platform: Siemens, Marconi, Philips, Toshiba, Hitachi to name a few.



The Secure Slide System includes:

- Secure Slide Attachment Extender Plate
- Secure Slide Wide Strap - 33" x 14"
- Secure Slide Narrow Strap - 33" x 6"
- Secure Slide Compatible Table Pad

Compatibility:

CT, MRI, Angiography, Nuclear Medicine

Applications:

Can be adapted to any imaging system

Product Features:

“Secure Slide Runner”



Benefits:

This unit provides the means for connecting the straps and other accessories to the table

“Secure Slide Connector”



Fixed to the straps and other attachments, this coupling allows the accessories to be slid up and down the runner

Options:

Standard sizes can be customized to exacting needs