


## **Siemens UROSKOP® Access Urology Drain Bag**

CFI Medical introduces a Urology Drain Bag designed specifically for the Siemens UROSKOP® Access. This unique design includes a vinyl flap that goes under the patient to direct fluids into the bag, a flexible 6' long hose and a rigid filter to allow for collection of tissue specimen.



<b><u>Part #</u></b>	<b><u>Description</u></b>	<b><u>Qty</u></b>	<b><u>List Price</u></b>
5442	Siemens UROSKOP® Access Drain Bag	Case of 20	\$225.00

**Options:**

CFI Medical manufactures Urology Drain bags to fit nearly every Urology system, including the OEC Uroview® 2000, 2600 & 2800, Liebel Flarsheim/Mallinckrodt Hydrovision ED, DR and HP and Dornier/Med Tec designs.

**Related Products:**

CFI Medical offers an entire line of Urology products including Urology Table Pads, Frames and Brackets, and a patented Urology Drain Collection Bag.