

MAXANT

MediPort Space-Saving Digital Viewing Work Station



Base unit with TS401 illuminator



Base unit with two 3MP monitors

Features:

- Designed and manufactured for digital medical applications
- Space-saving construction frees up counter space and takes up no more space than an illuminator
- Integrated digital imaging/analog illuminator systems option
- Choice of wall-mount or recess-mount
- Straightforward ease of operation, which includes a trackball (optional USB mouse)
- Configurations available include clinical review station or primary diagnostic work station
- Mounting hardware included for wall-mount system

Benefits:

- Direct wall replacement for MAXANT & other film viewers no additional wall space needed
- Medical grade: suited for use in isolated power systems such as in an OR, ER, ICU, CCU and other medical applications
- No change required in workflow or interaction with patients

Technical Support:

- Technical phone support 9 to 5 M-F E.S.T.
- Multiple tier technical support programs available
- Expanded support available 24/7

Warranty:

- One year parts and labor with base technical support.
- Extended warranty programs available on an annual basis.

Customer Service Philosophy:

We continue in the path that has built our business for more than three decades. We focus on superb service to our customers by continually exceeding our user's expectations.

TECHNICAL FEATURES

Base System:

- INTEL D915G series board with 2.93 gigahertz CELERON Processor
- Buss speed 533mhz
- L1 cache 16 kb
- L2 cache 256 kb
- 512 megabyte of RAM
- 40 gigabytes of hard drive (Seagate 7200.7 Barracuda)
- 1000T Ethernet
- CD//DVD ROM
- Keyboard with trackball
- 1.3 megapixel resolution portrait mode 19" color monitor with 450 NITS light output
- Monitor viewing angle 170 degrees vertical and horizontal
- Contrast ratio 600:1

Expanded systems:

- Base system RAM expanded to 2 gigabytes
- Base system hard drive expanded to 80 gigabytes
- 2 additional color 1.3 megapixel monitors
- 2 additional gray scale DICOM calibrated monitors either 2 or 3 megapixel resolution
- Combinations of MAXANT Tech Line Illuminators and digital imaging monitors up to 6 over 6 units.

Mechanical:

- Single digital imaging housing size is 14 x 21 x 5 ½ inches
- Single illuminator size is 14 x 21 x 3 ¼ inches
- All units mount in the same mounting configuration as the MAXANT Tech Line Illuminator product line, surface mount or recess mount

Electrical:

- 120 or 240 Volts at 50/60 cycles
- Medical grade line cord or units may be hard wired

Standards:

- UL60601-1, EN60601-1, NFPA 99 as a class 1 medical device
- FCC class B, CISPR 22
- CAN/CSA C22 2 NO.601.1

MAXANT

Contact information:

www.jzimaging.com

1-888-459-9729

JZ Imaging & Consulting, Inc. 1151 Westwood Dr. Willoughby, OH 44094