

Accommodates
BARIATRIC
Patients



BIODEX MEDICAL

NEW

Sound Pro Combination Table

Dual-purpose table accommodates ultrasound and echocardiography procedures.

- Ergonomically designed to help prevent musculoskeletal injuries
- Extra wide top with 500-lb capacity – perfect for bariatric patients
- Motor-controlled adjustments – including Fowler back, height, Trendelenburg motions and auto level
- Cardiac scanning cutout features a patented swing-out cushion that doubles as an armrest
- Adjustable sonographer cutout allows comfortable scanning positioning
- Vascular scanning arm board can be mounted to either side of the table
- Side rails pivot out of the way to allow unobstructed patient access
- Fully equipped with table accessories including hand and foot controllers, pivoting side rails, vascular scanning arm board, retractable stirrups, IV pole and paper dispenser



The Sound Pro Combination Table is a versatile table solution, combining the imaging features of both the Ultra Pro and Echo Pro, making it ideal for OB/GYN, general ultrasound and echocardiography. Fully equipped with ultrasound and echocardiography capabilities,

the Sound Pro is extra wide with 500-lb lift capacity and a vast height range down to 23", chosen for wheelchair transfers, patients with ambulatory difficulties and the bariatric population.

Ergonomically designed to suit both patient and sonographer, the table's height range allows seated or standing work positions while the pivoting side rails and adjustable sonographer cutout cushion allow comfortable shoulder, arm and hand positions when scanning. There are no barriers between sonographer and patient, allowing correct scanning posture to prevent musculoskeletal injury. The cardiac scanning cutout allows open access to the patient's left thorax area for an unobstructed apical approach. The motorized Fowler back is infinitely adjustable up to 80° via hand or foot controller. The table is so accommodating, stress-echo patients are off the treadmill, onto the table and into position for imaging quickly and comfortably.

For the patient, comfort is assured through lengthy procedures with cushioned arm supports, an amply thick mattress and a multitude of comfort-driven design details. For example, the unique swing-out cushion of the cardiac scanning cutout doubles as an arm support. The extra-wide tabletop comfortably accommodates the bariatric patient. In fact, the table is so accommodating, it reduces the total amount of time required to get a quality image.

The Sound Pro Combination Table is fully equipped... no additional options or accessories are required. Features include: 2"- and 3"-thick table padding covered in easy-to-clean Naugahyde®, two patient restraining straps with Velcro® hook and loop fasteners, five-inch swivel casters and a central locking system, accessible from either side of the table, to secure table position.

SPECIFICATIONS:

Dimensions: 71" l x 30" w (180.3 x 76.2 cm); 35.50" w (90.17 cm) with side rails

Height Adjustable: 23" to 39" (58.4 to 99.06 cm)

Tabletop: Three primary sections; torso section features sonographer and patient cutouts with Fowler back positioning; center section remains fixed and foot section drops down to 40° and 80° for stirrup access

Cardiac Scanning Cutout: with pivoting cushion that doubles as an arm board

Adjustable Sonographer Cutout: drops down or folds up for maximum access to the patient

Mattress: Torso section 2" (5 cm) thick, center and foot sections 3" (8 cm) thick; foam rubber pad covered with durable Naugahyde®

Upholstery Color: Graphite

Motions:

Height: 23" to 39" (58.4 to 99.06 cm)

Trendelenburg: 0° to ±15°

Fowler Back: 0° to 80° infinitely adjustable

Controls:

Hand Control: Activates height, Trendelenburg motions, Fowler back and auto level motions

Foot Control: Activates height, Trendelenburg motions, and Fowler back

Wheels: 5" (12.7 cm) swivel casters; central floor locking system

Patient Restraints: One body straps

Finish: Powder coat

Patient Lift Capacity: 500 lb (227 kg)

Weight: 350 lb (159 kg)

Power: 115 VAC or 230 VAC

Warranty: One year

Certifications: ETL and ETL® Listed to UL 2601 and CAN/CSA C22.2 No. 601.1-M90 Standards, EN 60601-1



058-710 Table, Sound Pro, Combination, 115 VAC\$7,995.00
Includes: hand and foot controllers, pivoting side rails, vascular scanning arm board, retractable stirrups, IV pole and paper dispenser.

058-715 Table, Sound Pro, Combination, 230 VAC
Includes: hand and foot controllers, pivoting side rails, vascular scanning arm board, retractable stirrups, IV pole and paper dispenser.

Replacement:

056-612 Table Paper, Hygenic, 12 rolls\$72.00
18" w x 225' l (46 cm x 69 m)

BIODEX