

H.D.R. Treatment Enclosure

TECHNICAL SPECIFICATIONS

	<i>Outpatient</i>	<i>Intra-Operative</i>	<i>Double Door</i>
<i>Dimensions</i>	10'6"w x 8'6"h x 7'0"d (75 square feet) (10' vertical height recommended)	10'6"w x 8'6"h x 7'0"d (75 square feet) (10' vertical height recommended)	10'6"w x 8'6"h x 7'0"d (75 square feet) (10' vertical height recommended)
<i>Internal Clearance</i>	8'9"w x 6'8"h x 5'6"d	8'9"w x 6'8"h x 5'6"d	8'9"w x 6'8"h x 5'6"d
<i>Manufacturing</i>	Shielded steel panels Lead batten strips	Shielded steel panels Lead batten strips	Shielded steel panels Lead batten strips
<i>Wall Panels</i>	3 1/4" total thickness with 2" of lead and 3/8" thick steel face plates	3 1/4" total thickness with 2" of lead and 3/8" thick steel face plates	3 1/4" total thickness with 2" of lead and 3/8" thick steel face plates
<i>Ceiling Panels</i>	3 3/8" total thickness with 2" of lead and 3/8" thick steel face plates	3 3/8" total thickness with 2" of lead and 3/8" thick steel face plates	3 3/8" total thickness with 2" of lead and 3/8" thick steel face plates
<i>Ports</i>	Physics (right front) Utilities (ceiling)	Physics (right front) Utilities (ceiling) Anesthesia (front)	Physics (right front) Utilities (ceiling)
<i>Interior Finish</i>	Metal stud interior Gypsum board Painted interior Suspended ceiling	Metal stud interior Operating room compliant	Metal stud interior Gypsum board Painted interior Suspended ceiling
<i>Lighting</i>	Overhead fluorescents Wall sconces on dimmers	Operating room compliant	Overhead fluorescents Wall sconces on dimmers
<i>Threshold</i>	1/8"th. Steel over 1/2"th. Pb	1/8"th. Steel over 1/2"th. Pb	1/8"th. Steel over 1/2"th. Pb
<i>Doors</i>	(1) equal pair 44"x84"x3"thick 2" thick lead with shielded astragal	(1) equal pair 44"x84"x3"thick 2" thick lead with shielded astragal	(2) equal pair 44"x84"x3"thick 2" thick lead with shielded astragal
<i>Floor</i>	Unshielded	Shielded	Unshielded
<i>Weight</i>	50,000 lbs.	61,000 lbs.	50,000 lbs.
<i>Compliance</i>	Meets all regulatory, safety requirements	Full operating room compliance	Meets all regulatory, safety requirements

Other details regarding manufacturing, materials, installation requirements, structural specifications, custom options, financing, physics, or regulatory, etc., can be obtained from your NELCO representative.

Lists of installed references can be available upon request.

04v1