

musicBox™

S O U N D S Y S T E M

Better Sound – Great Price



Features:

- Patented State-of-the-Art Magnetically Activated Stereo Transducer
- “Over-The-Ear” Headset for Maximum Noise Reduction
- Integrated Intercom
- CD Player
- Digital AM/FM Tuner
- Compatible with All MRI Systems, All Manufacturers and All Field Strengths

MEDRAD®



Better Sound – Great Price

Medrad's MusicBox System

Technologist control unit controls patient music volume and contains an integrated intercom with volume control and push to talk button for enhanced communication between technologist and the patient.

Stereo transducer converts electronic signals into high quality stereo sound for the patient. Our patented, magnetically activated transducer provides the highest quality stereo reproduction available.

Patient amplifier interfaces the stereo transducer to the audio source, controls all signal switching and provides power for the stereo transducer.

Deluxe headset provides full ear coverage and maximum noise attenuation.

Compact headset incorporates replaceable foam ear tips to provide noise attenuation and a low profile suitable for all imaging procedure and coils.

Installation kit provides cabling and filters necessary for professional system installation. Noise attenuation in excess of 100 dB at 100MHz ensures shielded room integrity and image quality.



Increased Efficiency

The *MusicBox* gives technologists an effective tool to increase throughput and reduce claustrophobic rejection.

Increased Image Quality

Music reduces patient anxiety and patient motion artifact, resulting in more effective use of scanner time and improved image quality.

Competitive Edge

With the *MusicBox* facility capabilities are increased improving the number of referrals and increasing patient satisfaction.

For More Information, Please Contact Your Local Medrad Representative, or Call 1-800-633-7231

www.medrad.com

info@medrad.com

MEDRAD®

Medrad, Inc. Corporate Office

One Medrad Drive, Indianola, PA 15051-0780 USA
Telephone: 412-767-2400 • Fax: 412-767-4128



Medrad is dedicated to continually improving the quality of all our products and services such that our customers' satisfaction, loyalty, and respect are unsurpassed. It is our policy to clearly understand and agree upon the valid requirements of the work we perform for our customers, both internal and external, and to pursue 100% conformance to those requirements...on time, every time.

Medrad, Inc. is ISO 9001/ISO 13485 certified.

Medrad, Quality for Life, Q and Design are federally registered trademarks, MusicBox is a trademark and Medrad is a service mark of Medrad, Inc. USA.

© 2003 Medrad, Inc. All Rights Reserved.

MB NA BA 202331 Rev. A