



# Burnett BiDirectional TMJ Device

Now acquire dynamic MR images of the TMJ as it opens and closes.



- Precisely opens and closes the jaws with simple one-handed operation
- Audible click indicates 1 mm movement of mouthpiece
- Greater bite area for reduced motion artifacts
- Improved strength and tracking method virtually eliminates device flex, for better images
- Scale indicates distance between mouthpieces, up to 55 mm
- Convenient disposable mouthpieces for single patient use

For more information call  
1-800-MEDRAD-1



**medRAD**<sup>®</sup>


# Burnett BiDirectional TMJ Device

Cat. No.	Description
TMJ 200	Burnett BiDirectional TMJ Device with three sets of disposable mouthpieces
TMJ 250	10 sets of disposable mouthpieces

## Warranty

Product is warranted to be free from defects in materials and workmanship for a period of twenty-four months from date of shipment. Please see the published warranty for details.

© 2001 Medrad, Inc. All Rights Reserved



**The Leading Image**

**MEDRAD®**

*The only source for Medrad Engineering, Medrad Quality, Medrad Service & Medrad Support.*

### Medrad, Inc. Corporate Office

One Medrad Drive  
Indianola, PA 15051-0780 USA  
**Telephone:** 1 (412) 767-2400  
**FAX:** 1 (412) 767-4128  
**OTHER:** 1 (800) 633-7231

### Medrad, Inc. Service Repair

One Medrad Drive  
Indianola, PA 15051-0780 USA  
**Telephone:** 1 (412) 767-2400  
**FAX:** 1 (412) 767-4126  
**OTHER:** 1 (800) 633-7237

### Medrad Europe B.V.

P.O. Box 205  
6190 AE Beek  
The Netherlands  
**Telephone:** +31 (0) 43-3585601  
**Fax:** +31 (0) 43-3656598

[www.medrad.com](http://www.medrad.com)

E-mail address for more information:  
[info@medrad.com](mailto:info@medrad.com)

Printed in the U.S.A.



Medrad is dedicated to continually improving the quality of all our products and services such that our customers' satisfaction, loyalty, and respect are unsurpassed. It is our policy to clearly understand and agree upon the valid requirements of the work we perform for our customers, both internal and external, and to pursue 100% conformance to those requirements...on time, every time.

Medrad reserves the right to modify the specifications and features described herein, or discontinue manufacture of the product described at any time without prior notice or obligation. Coils are compatible with all existing hardware and software at the time of release. Please contact your authorized Medrad representative for the most current information.

Medrad, Inc. is EN-ISO 9001/EN 46001 Certified  
Medrad Europe B.V. is EN-ISO 9002/EN 46002 Certified

Medrad, Quality for Life, Q and Design are federally registered trademarks. The Leading Image is a trademark, and Medrad is a service mark of Medrad, Inc. U.S.A.

88306-00-BA-01 Rev. A

